

Q300[®] M Mini



ORDER FORM

Billing Address Or Customer-ID	Shipping Address	Identical To Billing Address
C.-ID:	C.-ID:	
Name:	Name:	
Street:	Street:	
Town/PC:	Town/PC:	
E-Mail:	E-Mail:	
Phone:	Phone:	

Order Information			
Purchase Order No.:			
Order Date:	Day	Month	Year

Base Price **£6,250**

Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO LITE
Typical WC Weight	Without Batteries 80 kg	Seat Size Range	Width 385 - 510 mm, Depth 410 - 560 mm
	With Batteries 103 kg	Seat To Floor Height Range	410 mm - 490 mm

MWD Base Specifications

User Weight Capacity	136 kg	Drive Wheels	12", 13" or 14" Drive Wheels (Black)
Base Size	Width Min. 520 mm	Tyres / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 36 km	Battery Size	41 & 56 Ah Batteries
Drive Motors	2-Pole 6 /10 kph, 4-Pole 6/8 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

Q300 Drive Bases		£
<input type="checkbox"/>	QML000074 Q300 M Mini (MWD Base Chair)	6250

Standard Seat Size			£		
Seat Cushion Width		Seat Depth		JAY Seat Width - Only for JAY Cushions	
SW385	15" (385 mm) (1)	SD41	16" (410 mm) (3)	SW38	15" (38 mm)
SW41	16" (410 mm) (1)	SD435	17" (435 mm) (3)	SW40	15.75" (400 mm)
SW435	17" (435 mm) (1)	SD46	18" (460 mm) (3)	SW42	16.5" (420 mm)
SW46	18" (460 mm) (2)	SD485	19" (485 mm) (3)	SW44	17.3" (440 mm)
SW485	19" (485 mm) (2)	SD51	20" (510 mm) (3)	SW46	18" (460 mm)
SW51	20" (510 mm) (2)	SD535	21" (535 mm) (3) (4)	SW48	19" (480 mm)
		SD56	22" (560 mm) (3) (4)	SW50	19.7" (500 mm)

Backrest Frame Size		£
QML030154	Backrest, Width-Adjustable, 385 - 485 mm	0
QML030155	Backrest, Width-Adjustable, 460 - 510 mm	0

Seat Cushion Type		£
QML020098	Seat Cushion, Standard, Height 50 mm	0
QML020110	Omit Seat Cushion	0

Seat Rail		£
QML020123	Seat Rail, Short, For Seat Depth 410 - 485 mm	0
QML020124	Seat Rail, Long, For Seat Depth 510 - 560 mm	0

SEDEO Seat Cushion		£
QML020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	150
QML020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm - Comfort Cover	242

SEDEO Cushion Cover		£
QML020518	Seat Cushion Cover, Incontinence Protection	0
QML020080	Seat Cushion Cover, Comfort	0

SEDEO Seat Boards		£
QML020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	128

JAY Seat Cushions (6)		£
JAY020006	JAY Easy Visco With Flat Base	167
JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	592
JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	674

Seat Heights (1)		£
SH41	410 mm (2)	0
SH43	430 mm	
SH45	450 mm	0
SH47	470 mm	0
SH49	490 mm	0

(1) Comes with cushion width 410mm
(2) Comes with cushion width 485 mm
(3) Comes with cushion depth 50cm
(4) Not with SEDEO Lite back cushion QML03130-133

(5) Comes with seat cushion cover, incontinence protection
(6) Not all JAY Cushions are available in all sizes - for availability check our online configurator
(7) Only for cushions with seat depth 420 - 500 mm

Powered Seating Options		£
QML130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
QML130004	Posterior Tilt, Powered, 0° - 30° (4)	1051
QML130012	Posterior Tilt, Powered, -5° - 25° (4) (5)	1205
QML130007	Posterior Tilt, Powered, 0° - 50° (8) (9)	1218
QML130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (4)	1887

Backrest Frame Type		£
QML030147	Backrest, Tool Angle Adjustable, 86°, 90°, 94°, 98°, 102°, 106°	0
QML030148	Backrest, Manual Recline, Attendant Operated Only, 85° - 120°	303
QML030149	Backrest, Powered Recline, 85° - 120° (4)	545

Backrest Frame Option		£
QML030156	Backrest, Fixed	0
QML030098	Backrest, Flip Down, Quick Release	182

Backrest Cushion Type		£
QML030150	Sling Back, With Tension Adjustable Straps (Section A)	0

Sling Back Height		£
BH46	460 mm	0
BH485	485 mm	0
BH51	510 mm	0

JAY Shape Back		£
JAY080001	JAY Shape Back	237
JAY080006	JAY Shape Back, Oversized (7)	237

JAY Shape Backrest Height		£
BH47	470 mm	0
BH495	495 mm	0
BH52	520 mm	0
BH56	560 mm (10)	0
BH585	585 mm (10)	0
BH61	610 mm (10)	0

JAY Shape Back		£
JAY080010	JAY Shape Lateral Thoracic Pad, (LxH) 150 x 150 mm, Pair	237

SEDEO Lite Backrest		£
QML030130	SEDEO LITE, Back Cushion, Profiled, Height 500 mm, Comfort Cover	194
QML030131	SEDEO LITE, Back Cushion, Profiled, Adjustable Straps, Height 570 mm, Comfort Cover	278
QML030133	SEDEO LITE, Back Cushion, Deep Contour Active, Adjustable Straps, Height 570 mm, Comfort Cover	319

SEDEO Lateral Thoracic Supports		£
QML030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	157
QML030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	157

(1) Measured from front of the seat pan to the floor
(2) Not with powered tilt or lift & tilt
(3) Not with Seat lift & tilt QML130009
(4) Not with QML110036 Remote Joystick VR2. In case of R-Net QML110032 Actuator control & lights module, R-Net ISM-6L Required

(5) Not with seat height 415 or 440 mm
(7) Supplied 1" larger than seat width to allow greater immersion, not with seat width 485 / 510 mm
(8) Not With QML110036 Remote Joystick, VR2. In Case Of QML110004 Power Module, R-Net: QML110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

(9) Only With QML030147 Backrest, Tool Angle Adjustable
(10) Not with QML030148 Backrest, Manual Recline, Attendant Operated Only, 85° - 120° or QML030149 Backrest, Powered Recline, 85° - 120°

SEDEO Pro Backrest Size		£
Backrest Width		
BW38	380mm, For Effective Back Width 380 mm	0
BW42	420 mm, For Effective Back Width 380 mm - 420 mm	0
BW46	460 mm, For Effective Back Width 420 - 460 mm	0
BW50	500 mm, For Effective Back Width 460 - 500 mm	0
BW54	540 mm, For Effective Back Width 500 - 540 mm	0
BW56	560 mm, For Effective Back Width 520 - 560 mm	0
Backrest Height		
BH44	440 mm, For Effective Back Height 440 - 500 mm	0
BH50	500 mm, For Effective Back Height 500 - 560 mm	0
BH56	560 mm, For Effective Back Height 560 - 620 mm	0

SEDEO Pro Backrest Cushion		£
QML030116	SEDEO Pro, Back Cushion Flexible, Deep Contour, Active, Adjustable Straps	425
QML030117	SEDEO Pro, Back Cushion Flexible, Extra Deep Contour, Active, Adjustable Straps	446

SEDEO Pro Backrest Cushion Cover		£
QML030108	Back Cushion Cover, Comfort	0

Headrest Receiver		£
QML030120	Headrest Receiver, Backrest Frame Mount	0
QML030119	Headrest Receiver, Back Cushion Mount	0

Headrest Pad		£
QML030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	88
QML030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	103
QML030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	157
WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	151
WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	245

Headrest Stem (3)		£
QML030433	SEDEO, Headrest Mount	132
WHT030408	Whitmyer, Headrest Mount, LINX ²	165
WHT030409	Whitmyer, Headrest Mount, AXYS	136

JAY J3 Back Size		£
JAY030110	JAY J3 Back Size, 360 mm	0
JAY030111	JAY J3 Back Size, 410 mm	0
JAY030112	JAY J3 Back Size, 460 mm	0

JAY J3 Back Size, Mid Deep Contour		£
JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	793

JAY J3 Back Size, Deep Contour		£
JAY030021	JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	793

(1) Only use with SEDEO headrest mount
(2) Only use with Whitmyer headrest mount
(3) Choose a Whitmyer headrest mount and headrest receiver in case you would like to add a head control in the special controls section

JAY J3 Back, Spine Align Fit Kit		£
Spine Align Mini Fit Kit II		
JAY090032	JAY J3 Spine Align Mini Fit Kit II, 410 mm	59
JAY090033	JAY J3 Spine Align Mini Fit Kit II, 460 mm	59
JAY090034	JAY J3 Spine Align Mini Fit Kit II, 510 mm	59
Positioning Boxes		
JAY090035	JAY J3 Positioning Box For Back	183
JAY090036	JAY J3 Hardware Box For Back	600
Armrest Types		£
QML040039	Armrest, Flip Back	0
QML040123	Armrest, Reclining, Flip Back	224
QML040022	Armrest, Single Post, SEDEO LITE (1) (2)	121
QML040021	Armrest, Single Post, SEDEO PRO Removable (4)	
QML040019	Armrest, Single Post, SEDEO PRO Flip Back (3)	0
Side Guards		£
QML040033	Side Guard (7)	0
QML040034	Side Guard, Integral Foam (6) (9)	0
QML040035	Side Guard, Wedge Shaped, Padded (7) (9)	261
Arm Pads		£
QML040040	Arm Pad, Integral Foam, 380 mm, Right	0
QML040041	Arm Pad, Integral Foam, 380 mm, Left	0
QML040025	Arm Pad, Integral Foam, 400 mm, Right (4)	0
QML040026	Arm Pad, Integral Foam, 400 mm, Left (4)	0
QML040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	73
QML040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	73
QML040031	Arm Pad, Hemiplegia, Integral Foam, Right (5)	266
QML040032	Arm Pad, Hemiplegia, Integral Foam, Left (6)	266
Thigh Supports & Pommel (8)		£
QML020113	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right	182
QML020114	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left	182
QML020121	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right	182
QML020122	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left	182
QML030605	Abduction Pommel (9)	249
Centre Mount Legrest (Single Piece Footplate)		£
QML050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard	0
QML050087	Centre Calf Support, Integral Foam	0
QML050092	Calf Pad, Integral Foam, Pair	0
Power Centre Mount Legrest		£
QML050076	Centre Mount Legrest, Power 90° - 180°, 175 mm Extension	1438
Pad Options		
QML050072	Footplate, Dual Height and Angle Adjustment 260 x 150 mm	0
Pad Options		
QML050092	Calf Pad, Padded, Pair	0

SEDEO Lite Swing Away Legrest		£
QML050100	SEDEO LITE, Legrest, 70°, Swing Away	0
QML050105	SEDEO LITE, Legrest, Manual Elevating, Swing Away	351
Uni Swing Away Legrest		£
QML050115	Legrest, UNI, 80°, Swing Away	0
QML050142	UNI Amputee Supports, Padded, Right	133
QML050143	UNI Amputee Supports, Padded, Left	133
SEDEO Pro Swing Away Legrest		£
QML050015	SEDEO PRO, Legrest, 70°, Swing Away	0
QML050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	351
QML050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	653
Footplates		
QML050088	Footplate, For Seat Width 385 - 400 mm	0
QML050089	Footplate, For Seat Width 420 - 460 mm	0
QML050090	Footplate, For Seat Width 480 - 520 mm	0
QML050081	Footplates, Coupled (1)(2)	0
Hanger Receivers		
QML050170	SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thicknes 70mm And Higher	0
QML050171	SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thicknes 50mm And Lower	0
Calf Pad / Strap		
QML050082	Calf Support, Integral Foam, Depth-Adjustable, Pair	94
QML050083	Calf Support, Padded, Depth-Adjustable, Pair	128
QML050084	Heel Loop	0
QML050085	Calf Strap (2)	0
QML050200	Knee Pad, Integral Foam	67
VR2 Controls		£
QML110036	Remote Joystick, VR2	0
QML110037	Remote Joystick, VR2, 2 Actuators	0
QML110040	Remote Joystick, VR2, 2 Actuators And Lights	0
VR2 Power Module & Lighting Module		0
QML110001	Power Module, VR2, 70 A (4)	0
QML110002	Power Module, VR2, 90 A (5)	128
QML110030	Lighting Module, VR2 (6)	0
R-Net Controls		£
QML110015	Remote Joystick, R-Net, LED, Actuators	454
QML110016	Remote Joystick, R-Net, LED, Actuators And Lights	446
QML110017	Remote Joystick, R-Net, Colour LCD	587
QML110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	1174
QML110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	721
QML110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	908
R-Net Power Module & ISM		0
QML110003	Power Module, R-Net, 90 A	0
QML110004	Power Module, R-Net, 120 A	128
QML110032	Actuator Control & Lights Module, R-Net, ISM-6L	0
R-Net Gyro		0
QML110009	Traction Control System For R-Net (GYRO Technology)	394

(1) QML040033 Side Guard Included
(2) QML040040/41 Arm Pad, Integral Foam, 380 mm right/left included
(3) Not With QML040040-41 Arm Pad 380 mm, Integral Foam, Right Or Left; Not With Seat To Floor Height 415 mm
(4) Not With QML040022 Armrest, Single Post, SEDEO LITE

(5) Not With QML110052 Joystick Mount, Right or Left, QML090040 Tray Table, Swing Away, QML090046 Tray Table, Half-Side, Swing Away, QML090110 Rear View Mirror, Right

(6) Not With QML110053 Joystick Mount, Left, QML090040 Tray Table, Swing Away, QML090046 Tray Table, Half-Side, Swing Away, QML090118 Rear View Mirror, Left
(7) Only With QML040039, QML040123, QML040021 & QML040019
(8) Not with any centre mount legrest
(10) Only available with JAY Bracket

Joystick Mount		£
QML110050	Joystick Mounting, Fixed	0
QML110082	Joystick Mounting, Parallel, Swing Away	0
QML110056	Joystick Mounting, Parallel, Swing Away, Capsulated (1)	0

Joystick Mounting Location		£
QML110052	Joystick Mounting Location, Right	0
QML110053	Joystick Mounting Location, Left	0

Attendant Controls		£
QML110021	Dual Control	418
QML110020	Joystick Mounted As Attendant Control (2)	254

Attendant Mounting Location		£
QML110025	Attendant Mounting Location, Right	0
QML110026	Attendant Mounting Location, Left	0

Joystick Knobs		£
QML090505	Joystick Knob, Golf Ball, Ø 42 mm	39
QML090506	Joystick Knob, T-Bar, 82 mm Wide	39
QML090507	Joystick Knob, Soft Ball, Ø 70 mm	39
QML090508	Joystick Knob, Ball, Ø 35 mm	39
QML090045	Joystick Knob, Allround, Ball, Ø 48 mm	39
QML090509	Joystick Knob, Chin Cup, Ø 37 mm	39
QML090512	Joystick Knob, Stick, 80 x 15 mm	39
QML090511	Joystick Knob, Mushroom, Ø 45 mm	39

LED Lights & Indicators		£
QML110201	LED Light And Indicator StVZO Approved (3)	387

Motor Packages		£
QML180012	2-Pole Motor, 6 kph/ 4 mph, Performance	0
QML180015	2-Pole Motor, 10 kph/ 6.5 mph, Performance Package - Inc. LED Lights & Indicators, Rear View Mirror & Gyro (4)	473
QML180009	4-Pole Motor, 6 kph / 4 mph	303

C-Me - Drive at 4.8kph in full Lift (6)		£
QML110034	C-Me Function, Enabled (7)	0
QML110035	C-Me Function, Disabled (8)	0

Drive Wheels & Castor Wheels		£
QML070084	Drive Wheel, 12", Solid, Knobbly, Black	0
QML070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	61
QML070033	Drive Wheel, 14", Solid, Knobbly, Black	61
QML080213	Castor Wheel, 7", Solid, Black	0

(1) Only For QML050015 SEDEO PRO, Legrest, 70°, Swing Away
(2) Extra On A Set Of Swing Away Legrests, Attached Lose

(8) Not with any centre mount legrest
(9) Not With Seat To Floor Height 415 mm

Shroud Colours		£
QML100576	Red	0
QML100578	Blue	0
QML100579	White	0
QML100577	Black	0

Accent Colours		£
QML140005	Apple Green (RAL 110 70 77)	0
QML140004	Black (RAL 9005)	0
QML140006	Orange (RAL 060 60 80)	0
QML140007	Silver (RAL 9006)	0
QML140002	Red (RAL 040 40 67)	0
QML140003	Blue (RAL 260 40 45)	0
QML140001	White (RAL 9010)	0

Personalised Shrouds		£
QML140008	Grey Burst	237
QML140009	Red Burst	237
QML140010	Blue Burst	237
QML140011	Pure White Burst	237
QML140012	Flame Aesteroid Grey	237
QML140013	Flame Hot Red	237
QML140014	Flame Sharp Blue	237
QML140015	Flame Pure White	237
QML140016	Floral White	237
QML140017	Carbon Black	237

Batteries		£
QML160006	Battery, 41 Ah, Gel, Maintenance-Free	0
QML160013	Battery, 56 Ah, Gel, Maintenance-Free	121

Chargers		£
QML160100	Charger, Automatic, 220 V / 8 A, UK	0
QML160103	Charger, Automatic, 220 V / 10 A, UK	224

Lap Belt & JAY Belts		£
QML090128	Positioning Belt, Carseat Buckle, Short, 50 x 1700 mm	0
QML090123	Seat Rail Mounting Brackets For Positioning Belt	0
JAY090100	JAY Standard 2-Point Pelvic Belt, Side Release Buckle, Pad 190 x 40 mm	61
JAY090101	JAY Standard 2-Point Pelvic Belt, Side Release Buckle, Pad 240 x 40 mm	61
JAY090077	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 140 x 40 mm (1)	79
JAY090078	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 165 x 40 mm (1)	79
JAY090102	JAY V-Contour Pelvic Belt, Side Release Buckle, Pad Medium	109

Accessories		£
QML090110	Rear View Mirror, Right (6)	55
QML090118	Rear View Mirror, Left (6)	55
QML090009	Crutch Holder (7)	65

Tray Table (Mounted Opposite Joystick)		£
QML090124	Tray Table, Swing Away	327
QML090040	Tray Table, Width-Adjustable, Swing Away	327

(1) Not With QML040022 Armrest, Single Post

(2) Choose Also Remote Joystick

(3) QML110030 Lighting Module VR2 Or QML110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

(4) Only With 6 kph

(5) Highly Recommended With User Weight >100 kg

(6) If LED Light And Indicator Is Chosen, Then This Is Required

5-Switch Box		£
QML110043	#1 Ctrl+5 Switch Box With Free Assignable Functions	847
#1 Ctrl+5 Switch Box Mounting Location		
QML110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (1)	0
QML110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0
#1 Ctrl+5 Switch Type		
QML110175	#1 Ctrl+5 Switch Type With Push Buttons	0
QML110176	#1 Ctrl+5 Switch Type With Toggle Switches	0
Omni2 Input Control Module		£
QML110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared (2)	1145
QML110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared	1145
Omni 2 Module Mount		£
QML110119	OMNI2 Module, Fixed Mount, Right	0
QML110120	OMNI2 Module, Fixed Mount, Left	0
QML110123	OMNI2 Module, Swing Away Mount, Right	121
QML110124	OMNI2 Module, Swing Away Mount, Left	121
Hand Control		£
QML110301	Hand Control, Switch-It MicroGuide	2940
QML110302	Hand Control, Switch-It MicroPilot	3267
QML110304	Hand Control, Switch-It VersaGuide	1045
QML110305	Hand Control, Switch-It VersaGuide EZ (3)	1242
QML110306	Hand Control, Allround Joystick	1045
QML110307	Hand Control, Allround Joystick Light (4)	1242
QML110308	Hand Control, Allround Joystick Heavy Duty	1960
QML110309	Hand Control, Allround Joystick, R-Net Bus Enabled (5)	1111
QML110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled (5)	1307
QML110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (5)	2091
Hand Control Mounting		£
QML110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray (7)	653
Hand Control Mounting		
QML110350	Hand Control Mounting, Link-It To The Armrest, Fixed	194
QML110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	206
Hand Control Mounting Location		
QML110362	Hand Control Mounting Location, Right	0
QML110361	Hand Control Mounting Location, Left	0
Table Top Controls		£
QPL110143	Table Top Control, LCD, Centre Tray Mounted (4) (5) (6)	2091
Tray Table Mounting Location		
QPL110080	Tray Table Mounting Location, Left	0
QPL110081	Tray Table Mounting Location, Right	0

(1) Must also select R-Net Controls & Gyro

(2) (Choose Also QML090123 Seat Rail Mounting Brackets For Positioning Belt

(3) Only With QML1300009

(4) 4.5 kph Maximum Speed When Lifted (Tilt Angle <9°).

(5) 1.5 kph Maximum Speed When Lifted

Head Control (1)		£
Proportional Input Devices		
QML110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5619
QML110371	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	5619

Chin Control		£
QML110470	Chin Control, Switch-It MicroGuide	2940
QML110471	Chin Control, Switch-It MicroPilot	3267
QML110472	Chin Control, Switch-It VersaGuide	1045
QML110473	Chin Control, Switch-It VersaGuide EZ	1242
QML110474	Chin Control, Allround Joystick	1045
QML110475	Chin Control, Allround Joystick Light	1242
QML110476	Chin Control, Allround Joystick, R-Net Bus Enabled (2)	1111
QML110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (2)	1307

Chin Control Mounting		£
QML110480	Chin Control Mounting, Link-It To Chin Harness, Small (3)	589
QML110481	Chin Control Mounting, Link-It To Chin Harness, Large (3)	653
QML110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	523
QML110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	523
QML110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2614
QML110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2614
QML110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	458
QML110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	458

Link-It Extension Kit		£
QML110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	157
QML110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	196

Sip ,N Puff (8)		£
QML110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (9)	653

#1 Swiches - Section H (1)		£
QML110611	#1 Single Switch, Ø 64 mm, Black	113
QML110612	#1 Single Switch, Ø 64 mm, Red	113
QML110613	#1 Single Switch, Ø 64 mm, Green	113
QML110523	#1 Egg Switch – Black	124
QML110524	#1 MicroLite Switch – Black	183
QML110651	#1 Single Switch, Ø 30 mm, Black	113
QML110654	#1 Single Switch, Ø 30 mm, Blue	113
QML110653	#1 Single Switch, Ø 30 mm, Green	113
QML110652	#1 Single Switch, Ø 30 mm, Red	113
QML110655	#1 Single Switch, Ø 30 mm, Yellow	113
QML110534	#1 Twister Switch – Black	183
QML110535	#1 Twister Switch – Green	183
QML110536	#1 Twister Switch – Red	183
QML110537	#1 Satellite Switch – Black (7)	196
QML110538	#1 Satellite Switch – Green (7)	196
QML110539	#1 Satellite Switch – Red (7)	196

(1) Not Possible With QML110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)

(2) Display With Built In Infrared & Armrest Mount, Required For Specialty Input Devices
 (3) Not With QML090512 Joystick Knob, Stick, 80 x 15 mm, QML090511 Joystick Knob, Mushroom, Ø 45 mm, QML090506 Joystick Knob, T-Bar, 82 mm Wide

(4) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm

(5) OMNI Module Not Required

(6) Not In Combination With Same Entry Side QML090040 Tray Table, Swing Away, QML090046 Tray Table, Half-Side, Swing Away

(7) Not With QML030148 Backrest, Manual Recline, Attendant Operated Only, 85° - 120° & QML030149 Backrest, Powered Recline, 85° - 120°

#1 Switch Mounting		£
QML110560	#1 Switch Mounting, Link-It To The Hand Control (2)	196
QML110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (3)	196
QML110561	#1 Switch Mounting, Link-It To The Headrest	262
QML110563	#1 Switch Mounting, Link-It To The Seat Rail (4)	196
QML110565	#1 Switch Mounting, Link-It To Chin Boom (5)	196
QML110569	#1 Omit Switch Mounting (6)	0

#1 Switch Mount Location		£
QML110572	#1 Switch Mounting Location, Right	0
QML110571	#1 Switch Mounting Location, Left	0

#2 Swiches - Section H (1)		£
QML110631	#2 Single Switch, Ø 64 mm, Black	113
QML110632	#2 Single Switch, Ø 64 mm, Red	113
QML110633	#2 Single Switch, Ø 64 mm, Green	113
QML110088	#2 Egg Switch – Black	124
QML110089	#2 MicroLite Switch – Black	183
QML110671	#2 Single Switch, Ø 30 mm, Black	113
QML110674	#2 Single Switch, Ø 30 mm, Blue	113
QML110673	#2 Single Switch, Ø 30 mm, Green	113
QML110672	#2 Single Switch, Ø 30 mm, Red	113
QML110675	#2 Single Switch, Ø 30 mm, Yellow	113
QML110095	#2 Twister Switch – Black	183
QML110096	#2 Twister Switch – Green	183
QML110097	#2 Twister Switch – Red	183
QML110111	#2 Satellite Switch – Black (7)	196
QML110112	#2 Satellite Switch – Green (7)	196
QML110113	#2 Satellite Switch – Red (7)	196

#2 Switch Mounting		£
QML110098	#2 Switch Mounting, Link-It To The Hand Control (1)	196
QML110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (2)	196
QML110100	#2 Switch Mounting, Link-It To The Headrest	262
QML110104	#2 Switch Mounting, Link-It To The Seat Rail (3)	196
QML110105	#2 Switch Mounting, Link-It To Chin Boom (4)	196
QML110106	#2 Omit Switch Mounting (5)	0

#2 Switch Mount Location		£
QML110108	#2 Switch Mounting Location, Right	0
QML110107	#2 Switch Mounting Location, Left	0

Multi-Function Switches (7)		£
QML110550	Mini 2 Button (Stereo Plug) (8)	360
QML110551	Mini 2 Button (Dual Mono Plug) (8)	360
QML110552	Mini 4 Button InLine (DB9 Connector) (8)	604
QML110554	4-Way Toggle Switched Joystick	604

(1) For switch selection see section H. Choose a Whitmyer headrest mount.

(2) Omni 2 module required

(3) If you need switches mounted on the harness, you can only use All Round joystick

(4) Not with SW40, SW41, SW42, SW48, SW485

(5) Not with QPL040022 Armrest, Single Post, Sedeo Lite, QPL110309 Hand Control, Allround Joystick, QPL110310 Hand Control, Allround Joystick Light, QPL110311 Hand Control, Allround Joystick heavy Duty

(6) Only R-Net Controls, comes with joystick

(7) Only For Micro Series Joysticks

Accessories		£
QML110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	131
QML110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	169
QML110582	Extension Cable DB9 (1 m)	157
QML110583	Y-Cable, 2 Female Mono To 1 Male Stereo	196

Safety Devices		£
QML090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) (9)	915

Comments For Personalisations



Q300 M Mini with SEDEO LITE £

Base Price	
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date	Signature
------	-----------

Day Month Year

- (1) Not With HD Allround Joystick Or Space Disk Mount
- (2) Armrest Rail Mount
- (3) Knee, If 2 Switches Selected Only Left And Right Mount
- (4) No Left / Right Switch Selection Required
- (5) Switches Attached Lose

- (6) Only For Allround Joystick, Maximum 2 Selectable
- (7) Used In Combination With Assignable Button Programming
- (8) Attached Lose, When Ordered. No Link-It Fixation Possible
- (9) Not With VR2