

Breezy Ibis (Transit)

Invoice to:	
Name:	
Street:	
Town:	Post code:
E-Mail:	
Tel:	
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	Post code:
E-Mail:	
Tel:	Fax:
Marked for:	

Base Price
£2,695

Standard features - IBIS Pro Transit					
Y01117	Sedeo Pro Seat Frame	Y00135	12.5" solid rear wheels	Y01231	Side Supports PU
Y01126	Height adjustable push handles	Y01122	Folding backrest, tool operated	Y01255	Lap strap
Y00113	7" solid front wheels	Y01215	Removable armrests	Y00005	Anti-tip bracket, flip-up
Crash tested and Passed			Max User weight 160Kg		

UNIT-ID 0255 £

BASE FRAME			
SEDEO PRO SEAT FRAME			
Y01117	<input checked="" type="checkbox"/>	Sedeo pro seat frame	
Y01126	<input checked="" type="checkbox"/>	Height adjustable push handles	<i>(not with third party backrest)</i>
Y01124	<input type="checkbox"/>	Third Party Backrest - width adjustable	£458
Y01127	<input type="checkbox"/>	Angle adjustable push bar under 120kg	<i>(not included with third party/seat shell backrest frame/PowerSupport)</i> £289
SEAT HEIGHT (without cushion)			
Y00586	<input type="checkbox"/>	Low seat post for SH 32 - 42.5 / 12.5 - 16.5"	
Y00587	<input checked="" type="checkbox"/>	High seat post for SH 36.5 - 47cm / 14.5 - 18.5"	
Y00575	<input type="checkbox"/>	32cm measured to top of seat board	Y00581 <input type="checkbox"/> 41cm measured to top of seat board
Y00576	<input type="checkbox"/>	33.5cm measured to top of seat board	Y00582 <input type="checkbox"/> 42.5cm measured to top of seat board
Y00577	<input type="checkbox"/>	35cm measured to top of seat board	Y00583 <input type="checkbox"/> 44cm measured to top of seat board
Y00578	<input type="checkbox"/>	36.5cm measured to top of seat board	Y00584 <input type="checkbox"/> 45.5cm measured to top of seat board
Y00579	<input type="checkbox"/>	38cm measured to top of seat board	Y00183 <input type="checkbox"/> 47cm measured to top of seat board
Y00580	<input type="checkbox"/>	39.5cm measured to top of seat board	
Y00113	<input checked="" type="checkbox"/>	7" Solid front castors	
Y00085	<input checked="" type="checkbox"/>	Axle position middle 3	
Y00135	<input checked="" type="checkbox"/>	12.5" Solid rear wheels	
BRAKES & POWER SUPPORT			
Y00092	<input checked="" type="checkbox"/>	Footbrake	
Y00093	<input type="checkbox"/>	Footbrake for SW 54cm	
Y02047 (CSA)	<input type="checkbox"/>	Attendant brake for use 12.5" (not available with footbrake)	
Y00117	<input type="checkbox"/>	Power Support	£4,574
ADJUSTMENTS (tilt & recline, backrest options on next page)			
Y01512	<input checked="" type="checkbox"/>	Attendant operated, Gas spring tilt adjustable - max user weight 125kg. tilt 30°	
Y01513	<input type="checkbox"/>	Attendant operated, Gas spring tilt adjustable - Recommended for user weight 90kg - 160kg. tilt 30°	£107
Y00400	<input type="checkbox"/>	Adjustable tilt in space operated by user (60 - 90 kg) max. tilt 18° - RH	£473
Y00401	<input type="checkbox"/>	Adjustable tilt in space operated by user (90 - 125kg) max. tilt 18° - RH	£473
Y02039 (CSA)	<input type="checkbox"/>	Negative Tilt. tilt 5°	£153
Y01119	<input type="checkbox"/>	Attendant recline. recline 89-128°	£289
Y01121	<input type="checkbox"/>	Manually reclining backrest, operated by the user - RH	£289
Y01120	<input type="checkbox"/>	Manually reclining backrest, operated by the user - LH	£289

				£
SEAT SIZE				(adjustability)
Y01139	<input type="radio"/>	42 x 50cm (16" x 20")	16" - 21" Seat Width x 17" - 20" Depth	
Y02644	<input type="radio"/>	42 x 56 cm (16" x 22")	16" - 21" Seat Width x 19" - 22" Depth	
Y01142	<input type="radio"/>	46 x 50cm (18" x 20")	18" - 21" Seat Width x 17" - 20" Depth	
Y02645	<input type="radio"/>	46 x 56 cm (18" x 22")	18" - 21" Seat Width x 19" - 22" Depth	
Y01145	<input type="radio"/>	50 x 50cm (20" x 20")	19" - 21" Seat Width x 17" - 20" Depth	
Y02646	<input type="radio"/>	50 x 56 cm (20" x 22")	19" - 21" Seat Width x 19" - 22" Depth	
Y01148	<input type="radio"/>	54 x 50 cm (21" x 20")	21" Seat Width x 17" - 20" Depth	
Y02647	<input type="radio"/>	54 x 56 cm (21" x 22")	21" Seat Width x 19" - 22" Depth	
SEAT CUSHIONS				
Y01159	<input checked="" type="checkbox"/>	Board Seat Pan seat with front (50 mm) & depth adj. rear lip		
Y01156	<input type="checkbox"/>	Stabilizer seat, 90mm thick		£153
Y01166	<input type="checkbox"/>	Abduction support (pommel)	<i>(not included with central legrest)</i>	£229
BACKRESTS				
BACKREST SIZE				
Y01173	<input type="radio"/>	42cm Back width - Standard (set at 50cm, adjustable to 58cm)		
Y01174	<input type="radio"/>	42cm Back width - Tall (set at 56cm, adjustable to 64cm)		
Y01176	<input type="radio"/>	46cm Back width - Standard (set at 50cm, adjustable to 58cm)		
Y01177	<input type="radio"/>	46cm Back width - Tall (set at 56cm, adjustable to 64cm)		
Y01179	<input type="radio"/>	50cm Back width - Standard (set at 50cm, adjustable to 58cm)		
Y01180	<input type="radio"/>	50cm Back width - Tall (set at 56cm, adjustable to 64cm)		
Y01182	<input type="radio"/>	54cm Back width - Standard (set at 50cm, adjustable to 58cm)		
Y01183	<input type="radio"/>	54cm Back width - Tall (set at 56cm, adjustable to 64cm)		
BACKREST OPTIONS				
Y01122	<input checked="" type="checkbox"/>	Folding backrest, tool operated	<i>(Standard option if nothing selected)</i>	
Y01123	<input type="checkbox"/>	Folding backrest, thumbwheel		£92
BACKREST CUSHIONS				
Y01190	<input type="radio"/>	Flexi Deep Contour, Passive, Adj. Straps	<i>(not compatible with thoracic supports)</i>	
Y01191	<input type="radio"/>	Flexi Deep Contour, Active, Adj. Straps	<i>(not compatible with thoracic supports)</i>	
Y01520	<input type="radio"/>	Flexi Extra Deep Contour, Active, Adj. Straps	<i>(not compatible with thoracic supports)</i>	
Y01185	<input type="radio"/>	Profiled Backrest	<i>(compatible with thoracic supports)</i>	
Y01193	<input checked="" type="checkbox"/>	Removable incontinence backrest cushion cover (included with all backrests)		
Y01194	<input type="checkbox"/>	Comfort cover (Breathable cover including incontinence inner liner)		£122
Y01195	<input type="checkbox"/>	Thoracic support (not included with third party or seat shell backrest frame)	<i>(Use with padded backrest)</i>	£442
LEGRESTS				
Y01198	<input type="radio"/>	Swingaway legrests		
Y01202	<input checked="" type="checkbox"/>	Calf Strap		
Y01204	<input type="checkbox"/>	Depth adj. calf support, PU, set		£198
Y01209	<input type="radio"/>	Centre Mount Legrest		
Y01210	<input type="radio"/>	Centre Mount Articulating Legrest	<i>(from seat height 35cm)</i>	
Y01212	<input checked="" type="checkbox"/>	Central calf support		
Y01213	<input type="checkbox"/>	Central calf support with side support		£92
Y01199	<input type="checkbox"/>	Manually elevating legrests	<i>Price on order will be made up of £236 on legrest and £180 on support pad</i>	£458
Y01204	<input checked="" type="checkbox"/>	Depth adj. calf support, PU, set		£198
Y01205	<input checked="" type="checkbox"/>	Knee pads (not available with padded side supports, wedge-shape)		
AMPUTEE SUPPORTS				
Y01208	<input type="checkbox"/>	Amputee support, right	<i>(not available with centre mount articulating legrest)</i>	£183
Y01207	<input type="checkbox"/>	Amputee support, left	<i>(not available with centre mount articulating legrest)</i>	£183

ARMRESTS & ARMPADS			
Y01216	●	Removable armrests	
Y01223	○	RIGHT HAND - Armrest pad 30cm	
Y01222	○	LEFT HAND - Armrest pad 30cm	
Y01221	○	RIGHT HAND - Armrest pad 40cm	
Y01220	○	LEFT HAND - Armrest pad 40cm	
Y01548	○	RIGHT HAND - Armrest pad wide 45cm <i>(not i.c.w. SW 54 cm, not i.c.w. 24" rear wheels)</i>	
Y01547	○	LEFT HAND - Armrest pad wide 45cm <i>(not i.c.w. SW 54 cm, not i.c.w. 24" rear wheels)</i>	
Y01504	○	Raised Armrests +10 cm (only i.c.w. removable armrests with quick-lock)	
Y01227	□	RIGHT HAND - Hemiplegia armrest and pad	£244
Y01226	□	LEFT HAND - Hemiplegia armrest and pad	£244
SIDE GUARDS			
Y01231	●	Side Supports PU (SW -2cm) <i>(Standard Option if nothing selected)</i>	
Y01232	□	Padded side supports, wedge-shape (SW -8cm)	£289
Y01234	□	Thigh supports <i>(available with Standard Side support but not available with padded wedge side support)</i>	£289
HEADRESTS (please pick headrest & mounting)			
Y01235	○	Headrest with PU Cover	
Y01236	□	Headrest flexible	£107
Headrest Mounting			
Y01239	○	Headrest socket on backrest frame <i>(Standard Option if nothing selected)</i>	
Y00333	○	Headrest socket on backrest frame - only with third party backrest	
Y01240	□	Width adjustable headrest bracket <i>(allows greater range of lateral adjustment)</i>	£107
ACCESSORIES			
Y00005	●	Anti-tip bracket, flip up	
Y01255	●	Lap strap	
Y00012	□	Cane/crutch holder for 2 canes	£67
WORKTOPS			
Y01249	□	Slide on Tray <i>(Fixed acrylic with two insert tubes, abdominal recess 11cm)</i>	£289
Base Price			£2,695
Option price			£
Total price			£

 <p>Armrest assembly (QR removable)</p>	 <p>Armpad 30cm / 40cm</p>	 <p>Hemiplegic armpad</p>
 <p>Side Supports PU</p>	 <p>Padded side supports, wedge-shape (SW -8cm)</p>	 <p>Thigh support</p>
 <p>Passive backrest - 14cm contour</p>	 <p>Active backrest - 14cm contour</p>	 <p>Active backrest - 21cm contour</p>
 <p>Thoracic Supports swing-away</p>	 <p>Abduction support (pommel) removable</p>	 <p>Amputee support</p>
 <p>Headrest PU (Std)</p>	 <p>Flexible Headrest</p>	 <p>Width adjustable headrest mount</p>
 <p>Swing away legrests</p>	 <p>Centre Mount - manually articulating and elevating legrest</p>	 <p>Manual ELR</p>
 <p>Central calf support with side support</p>	 <p>Third party Backrest frame</p>	 <p>Upholstery: Standard / Breathable</p>