



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

MARCH 2023
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Küschall Compact 2.0

Standard Features

- Max user weight 130kg/20.5 stone
- Crash tested ISO 7176-part 19
- Hydroformed aluminium, active, folding wheelchair, with swingaway or fixed frame
- Multiple adjustments including front and rear seat heights and centre of gravity
- Hydroformed swingaway legrests, choice of 70° and 80° angles
- Hydroformed fixed frame, choice of 70° or 80° angles
- Height adjustable flip up plastic footrests
- Choice of seat widths 280mm - 500mm (11" - 20")
- Choice of 10 seat depths 320mm - 500mm (12" - 20")
- Height adjustable fixed angle backrest set at 90°
- Standard light black nylon upholstery, wipeable
- Standard push handles
- Fixed carbon clothesguards
- 5" comfort solid front castors with black tyres
- Küschall starec q/r rear wheels with light pneumatic tyres and aluminium handrims
- Pull to lock brakes
- Choice of frame colours



Standard Base price: £1,920

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DCA0012	KÜSCHALL COMPACT 2.0	£1,920
DCA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable	POA

FRAME*			
DCA1100		Swing away Front frame	(NC 1900, 1902, 3040, 3140) ○
DCA1110		Fixed front frame 70°	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130) ○
DCA1120		Fixed front frame 80° (Minimum seat angle, front and rear seat heights is 40 mm with 5" front castors)	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650) ○
DCA1130		Fixed front frame tapered 70° Additional 10 day leadtime	(NC 1300, 1305, 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130) ○
DCA1140		Fixed front frame tapered 80° (minimum seat angle, front and rear seat heights is 40 mm with 5" front castors) Additional 10 day leadtime	(NC 1300, 1305, 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 2748, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650) ○
SEAT WIDTH*			
DCA1300		280 mm Single cross brace	(NC 2210, 2490, 2710, 2720, 2748, 2750, 2748, 2750, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○
DCA1305		300 mm Single cross brace	(NC 2210, 2490, 2710, 2720, 2748, 2750, 2748, 2750, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○
DCA1310		320 mm Single cross brace	(NC 2490, 2710, 2720, 2748, 2750, 2748, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○
DCA1315		340 mm Single cross brace	(NC 2490, 2710, 2720, 2748, 2750, 2748, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○
DCA1320		360 mm Single cross brace	(NC 2490, 2710, 2748, 3020, 3030, 3120, 3130) ○
DCA1325		380 mm Single cross brace	(NC 2490, 3020, 3030, 3120, 3130) ○
DCA1330		400 mm Single cross brace	◆
DCA1335		420 mm Single cross brace	○
DCA1340		440 mm Double cross brace	○
DCA1345		460 mm Double cross brace	○
DCA1350		480 mm Double cross brace	○
DCA1355		500 mm Double cross brace	(NC 2210) ○
SEAT DEPTH*			
DCA1400		320 mm	(NC 2155, 2160, 2710, 2720, 2748, 2750, 5650) ○
DCA1405		340 mm	(NC 2155, 2160, 2710, 2720, 2748, 5650) ○
DCA1415		360 mm	(NC 2710, 2720, 2748, 5650) ○
DCA1425		380 mm	(NC 2720, 5650) ○
DCA1435		400 mm	◆
DCA1445		420 mm	○
DCA1460		440 mm	○
DCA1470		460 mm	○
DCA1480		480 mm	(NC 2748) ○
DCA1490		500 mm	(NC 2710, 2748) ○

BACK HEIGHT Transportation kit is only possible with backrest heights of 405mm and above*

DCA1615		300 mm	(NC 2140, 2210, 2595, 7600)	○
DCA1620		315 mm	(NC 2140, 2210, 2595, 7600)	○
DCA1625		330 mm	(NC 2140, 2210, 2595, 7600)	○
DCA1630		345 mm	(NC 2210, 2595, 7600)	○
DCA1635		360 mm	(NC 2595, 7600)	○
DCA1640		375 mm	(NC 2595, 7600)	○
DCA1645		390 mm	(NC 2595, 7600)	○
DCA1650		405 mm		○
DCA1655		420 mm		◆
DCA1660		435 mm		○
DCA1665		450 mm		○
DCA1670		465 mm		○
DCA1675		480 mm		○
DCA1680		495 mm		○
DCA1685		510 mm		○

FRONT SEAT HEIGHT. Seat heights 380mm TO 430mm offers a unique front frame design to achieve a lower front seat height, only available for the Compact Swing option*

DCA1700		370 mm	(NC 1810, 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1932, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1705		380 mm	(NC 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1710		390 mm	(NC 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1715		400 mm	(NC 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1720		410 mm	(NC 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130)	○
DCA1725		420 mm	(NC 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1730		430 mm	(NC 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1735		440 mm	(NC 1845, 1850, 1855, 1860, 1865, 1870, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1740		450 mm	(NC 1850, 1855, 1860, 1865, 1870, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1745		460 mm	(NC 1855, 1860, 1865, 1870, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1750		470 mm	(NC 1860, 1865, 1870, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1755		480 mm	(NC 1805, 1865, 1870, 1954, 1956, 1958, 1960, 1962)	○
DCA1760		490 mm	(NC 1805, 1810, 1956, 1958, 1960, 1962)	○
DCA1765		500 mm	(NC 1805, 1810, 1815, 1958, 1960, 1962)	◆
DCA1770		510 mm	(NC 1805, 1810, 1815, 1820, 1960, 1962)	○
DCA1775		520 mm	(NC 1805, 1810, 1815, 1820, 1825, 1962)	○
DCA1780		530 mm	(NC 1805, 1810, 1815, 1820, 1825, 1830)	○

REAR SEAT HEIGHT*				
DCA1805		370 mm	(NC 6000, 6010, 6410, 6420, 6430)	○
DCA1810		380 mm	(NC 6000, 6010, 6410, 6420, 6430)	○
DCA1815		390 mm	(NC 6420, 6430)	○
DCA1820		400 mm	(NC 6420, 6430)	○
DCA1825		410 mm	(NC 6420, 6430)	○
DCA1830		420 mm		○
DCA1835		430 mm		○
DCA1840		440 mm		○
DCA1845		450 mm		◆
DCA1850		460 mm		○
DCA1855		470 mm		○
DCA1860		480 mm	(NC 6400)	○
DCA1865		490 mm	(NC 6400)	○
DCA1870		500 mm	(NC 6000, 6010, 6400, 6410)	○
LOWER LEG LENGTH. Cushion height must be deducted*				
DCA1900		200 mm	(NC 3030, 3130)	○
DCA1902		210 mm	(NC 3030, 3130)	○
DCA1904		220 mm	(NC 3030, 3130)	○
DCA1906		230 mm	(NC 3030, 3130)	○
DCA1908		240 mm	(NC 3030, 3130)	○
DCA1910		250 mm	(NC 3030, 3130)	○
DCA1912		260 mm	(NC 3030, 3130)	○
DCA1914		270 mm	(NC 3030, 3130)	○
DCA1916		280 mm	(NC 3030, 3130)	○
DCA1918		290 mm	(NC 3030, 3130)	○
DCA1920		300 mm	(NC 3030, 3130)	○
DCA1922		310 mm	(NC 3030, 3130)	○
DCA1924		320 mm	(NC 3030, 3130)	○
DCA1926		330 mm	(NC 3030, 3130)	○
DCA1928		340 mm	(NC 3030, 3130)	○
DCA1930		350 mm	(NC 3030, 3130)	○
DCA1932		360 mm	(NC 3030, 3130)	○
DCA1934		370 mm	(NC 3030, 3130)	○
DCA1936		380 mm		○
DCA1938		390 mm		○
DCA1940		400 mm		○
DCA1942		410 mm		○
DCA1944		420 mm		○
DCA1946		430 mm		◆
DCA1948		440 mm		○
DCA1950		450 mm		○
DCA1952		460 mm		○
DCA1954		470 mm		○
DCA1956		480 mm		○
DCA1958		490 mm		○
DCA1960		500 mm		○
DCA1962		510 mm	(NC 3000, 3100, 3010, 3110)	○

BACKREST TYPE*			
DCA2110		Fixed straight backrest tubes set at 90°	(NC 2155, 2160, 2165, 2170, 2175, 2180) ◆
DCA2120		Fixed backrest with lumbal bend angled at 7°	(NC 2155, 2160, 2165, 2170, 2175, 2180) ○
DCA2130		Angle adjustable back (angle must be selected from the options below)	(NC 0000) £192
DCA2140		Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below)	(NC 0000, 2180, 2430, 3520, 3620) £196
BACKREST ANGLE*			
DCA2155		Angle adjustable backrest, adjusted 82°	○
DCA2160		Angle adjustable backrest, adjusted 86°	○
DCA2165		Angle adjustable backrest, adjusted 90°	○
DCA2170		Angle adjustable backrest, adjusted 94°	○
DCA2175		Angle adjustable backrest, adjusted 98°	○
DCA2180		Angle adjustable backrest, adjusted 102°	○
DCA0000		Without angle adjustment. select with options 2110 fixed and 2120 lumbal backrests, supplied with fixed angle	○
BACKREST UPHOLSTERY*			
DCA2200		Black Standard Light Upholstery, without Tension Adjustment	(NC 2400, 2430, 2595) ◆
DCA2210		Lightweight back upholstery, no tension adjustment	(NC 2430, 2595) £66
DCA2220		Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener	£130
DCA2230		Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener	£212
DCA2240		Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required	(NC 2430, 2595) ○
PUSH HANDLES*			
DCA2400		Without Push Handles	(NC 2490) ○
DCA2410		Push Handles Fixed Short/Mini Grip	(NC 2490, 6580) ○
DCA2420		Push Handles Fixed Long Grip	◆
DCA2430		Push handles height adjustable, rear set (maximum height of push handles is rear seat height plus 600mm)	(NC 2490, 3520, 3620) £241
DCA2440		Push Handles Height Adjustable, Integrated in backrest	(NC 2490, 2595) £208
DCA2450		Push Handles Foldable	(NC 2490, 2595, 6580) £188
HEADREST			
DCA2490		Headrest	(NC 2595, 6580) £160
BACKREST STABILIZATION BAR			
DCA2595		Backrest Stabilization Bar, Swingaway and Removable	£156
SEAT UPHOLSTERY*			
DCA2610		Standard seat upholstery with 1 central band (loop)	○
DCA2620		Seat upholstery with 1 central (loop) and 2 side bands (hook)	◆
SEAT CUSHION			
DCA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850) £63
DCA2710		Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850) £115
DCA2720		Matrix Libra Cushion, Black Smooth Cover 70mm High - Very High Risk	£390
DCA2748		Contour Visco NG Cushion 100mm High - High Risk	(NC 2810, 2820, 2830, 2840, 2850) £161
DCA2750		Vicair Active O2 lightweight cushion 60mm High - High risk	(NC 2810, 2820, 2830, 2840, 2850) £495

SEAT BAG			
DCA2800		Seat Bag, Click in and removable	£70
MATRX LIBRA CUSHION ACCESSORIES			
DCA2810		Fluid supplement sac (Pair)	£42
DCA2820		Tapered Wedge Pair	£29
DCA2830		Abductor	£22
DCA2840		Lateral Side Wedge	£19
DCA2850		Anterior Front Wedge	£26
LEGREST TYPE RH*			
DCA3000		Legrest 70° RH (NC 3110, 9999)	◆
DCA3010		Legrest 80° RH (NC 3100, 5650, 9999)	○
DCA3020		Amputee legrest RH (NC 3280, 3300, 3305)	£131
DCA3030		Angle adjustable/Elevating legrest - RH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305)	£113
DCA3040		Without Legrests RH, Select with fixed front frames	○
LEGREST TYPE LH*			
DCA3100		Legrest 70° LH (NC 9999)	◆
DCA3110		Legrest 80° LH (NC 5650, 9999)	○
DCA3120		Amputee legrest LH (NC 3280, 3300, 3305)	£131
DCA3130		Angle adjustable/Elevating legrest - LH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305)	£113
DCA3140		Without Legrests LH, Select with fixed front frames	○
FOOTREST*			
DCA3260		2 pc. Plastic flip up	◆
DCA3270		2 pc. Carbon angle and depth adjustable (NC 3310)	£316
DCA3280		1 pc. Carbon flip-up to the left, angle and depth adjustable (NC 3310)	£159
DCA9999		No footrest (Select when ordering Amputee brackets both sides)	○
CALF STRAP			
DCA3300		Calf strap mounted on the frame	○
DCA3305		Calf strap mounted on bracket	○
HEEL LOOPS			
DCA3310		Heel loops	○
CLOTHESGUARDS/MUDGUARDS*			
DCA3405		Composite clothes guards, removable (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450)	◆
DCA3410		Carbon clothes guards, fixed (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450)	○
DCA3430		Carbon mudguards, fixed (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 7600, 8450)	£235
DCA3440		Carbon mudguards, removable (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6140, 6150, 8450)	£292
DCA6666		Without clothesguards or mudguards, select only with options 3500, 3510, 3530, 3540, 3600, 3610, 3630 or 3640 Armrests supplied with clothesguards (NC 3520, 3620)	○

ARMREST RH*			
DCA3500		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH	(NC 6150) ○
DCA3510		Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH	(NC 6150) ○
DCA3520		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included	(NC 8450) £113
DCA3530		Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH	(NC 6130, 6140, 6150) £201
DCA3540		Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH	(NC 6130, 6140, 6150, 8450) £98
DCA7777		Without armrest, select with Clothesguard/Mudguard options	(NC 3500, 3510, 3530, 3540) ○
ARMREST LH*			
DCA3600		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH	○
DCA3610		Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH	○
DCA3620		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH	(NC 8450) £113
DCA3630		Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH	(NC 6130, 6140, 6150) £201
DCA3640		Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH	(NC 6130, 6140, 6150, 8450) £98
DCA8888		Without armrest, select with Clothesguard/Mudguard options	(NC 3600, 3610, 3630, 3640) ○
FRAME COLOUR*			
DCA4000		Space Black (RAL 9005)	○
DCA4010		Anthracite Grey (RAL 7016)	○
DCA4020		Pure White (RAL 9016)	(NC 4220) ○
DCA4030		Carmine Red (RAL 3002)	○
DCA4040		Azure Blue (RAL 5017)	○
DCA4060		Flashy Yellow (RAL 1016) (Additional 10 day lead time)	£29
DCA4070		Atomic Violet (Additional 10 day lead time)	£29
DCA4100		Black Glossy (RAL 9005) (Additional 10 day lead time)	£29
DCA4130		Silver Grey (RAL 9006) (Additional 10 day lead time)	£29
DCA4140		Royal Blue (Additional 10 day lead time)	£29
DCA4150		Ruby Red (RAL 3003) (Additional 10 day lead time)	£29
DCA4170		Dynamite Orange (Additional 10 day lead time)	£29
DCA4180		Forest Green (RAL 6029) (Additional 10 day lead time)	£29
DCA4050		Personal frame colour. Provide standard RAL colour code of your choice	£295
FRAME DECAL TYPE*			
DCA4190		Bold Decal (Full coloured transfer on frame)	(NC 4300, 4310, 1111) £29
DCA4195		Tag (Small Küschall brand transfer on front of frame)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) ◆
DCA4197		Without Tag (Kuschall branding only)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) ○

BOLD DECAL COLOUR*			
DCA4200	<input type="checkbox"/>	Black	○
DCA4210	<input type="checkbox"/>	Dark Grey	○
DCA4220	<input type="checkbox"/>	Light Grey	○
DCA4230	<input type="checkbox"/>	Red	○
DCA4240	<input type="checkbox"/>	Blue	○
DCA4250	<input type="checkbox"/>	Gold	○
DCA1111	<input type="checkbox"/>	No bold colour, select only with Option DCA4195 Tag (NC 4190) or DCA4197 Without Tag	○
TAG COLOUR*			
DCA4300	<input type="checkbox"/>	Black	○
DCA4310	<input type="checkbox"/>	Red	○
DCA2222	<input type="checkbox"/>	No tag colour, select only with Option DCA4190 Bold (NC 4195) Decal or DCA4197 Without Tag	○
FRONT CASTORS*			
DCA5050	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	○
DCA5100	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	○
DCA5150	<input type="checkbox"/>	100mm x 25mm (4" x 1") Starec Black hub and tyre	£79
DCA5200	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£88
DCA5250	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£88
DCA5300	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	○
DCA5350	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	◆
DCA5400	<input type="checkbox"/>	125mm x 25mm (5" x 1") Starec Black hub and tyre	£79
DCA5450	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£88
DCA5500	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£88
DCA5550	<input type="checkbox"/>	150mm x 30mm (6" x 1.2") Comfort Black hub and tyre	○
DCA5650	<input type="checkbox"/>	175mm x 30mm (7" x 1.2") Comfort Black hub and tyre	○
ONE ARM DRIVE			
DCA6000	<input type="checkbox"/>	One Arm Drive system Right side (with light wheel 7 bar tyres, alu handrims)	(NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6600, 6630, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) £458
DCA6010	<input type="checkbox"/>	One Arm Drive system Left side (with light wheel 7 bar tyres, alu handrims)	(NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6600, 6630, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) £458
REAR WHEEL FIXATION*			
DCA6020	<input type="checkbox"/>	Standard fixation, 4 rear wheel positions	(NC 6140, 6150) ◆
DCA6040	<input type="checkbox"/>	Rear wheel extension, 6 rear wheel positions	(NC 6100, 6580) £63
REAR WHEEL POSITION. Position 1 is most active and position 6 is most stable*			
DCA6100	<input type="checkbox"/>	1st Position	(NC 6580, 6840) ○
DCA6110	<input type="checkbox"/>	2nd Position	○
DCA6120	<input type="checkbox"/>	3rd Position	◆
DCA6130	<input type="checkbox"/>	4th Position	(NC 6580) ○
DCA6140	<input type="checkbox"/>	5th Position	○
DCA6150	<input type="checkbox"/>	6th Position	○
QUAD QUICK RELEASE AXLE			
DCA6300	<input type="checkbox"/>	Quad quick release/tetra clip	(NC 6550, 6580) £81

REAR WHEEL SIZE*			
DCA6400		22"	(NC 6520, 6550, 6600, 6630, 6660, 6760, 6710, 6720, 6740, 6750, 6760, 6770, 4444) ○
DCA6410		24"	◆
DCA6420		25"	(NC 6510, 6580, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900, 4444) ○
DCA6430		26"	(NC 6520, 6550, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6840, 4444) ○
REAR WHEEL TYPE*			
DCA6510		22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim	◆
DCA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6660) £173
DCA6550		24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black	(NC 6660) £580
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall	(NC 6000, 6010, 6660) £462
SWO0909		24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall	(NC 6000, 6010, 6660) £698
DCA6580		22", 24", 26" Wheels with Attendant Drum Brake	(NC 8680) £348
DCA3333		select with options 6000 and 6010 One Arm Drive wheels, specific wheel supplied	○
REAR WHEEL TYRE*			
DCA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DCA6625		24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.	£95
DCA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	£92
DCA6660		24" x 1 3/8" Solid Profile Black Tyre	(NC 7230) £63
DCA4444		Select with options 6000 and 6010 One Arm Drive wheels, supplied with specific tyre	○
HANDRIM*			
DCA6700		22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	£212
DCA6710		24", 26" Steel, short tab with 1 hole, mounted at 3cm distance	£59
DCA6720		24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	£182
DCA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip	£304
DCA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands	£304
DCA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion	£270
DCA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim	£338
DCA6780		22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
DCA6790		22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	£40
DCA5555		select with options 6000 and 6010 one arm drive wheels, supplied with specific handrims	○

CAMBER ANGLE*			
DCA6820		Camber 0°	◆
DCA6840		Camber 3°	(NC 7230) ○
SPOKEGUARD			
DCA6900		22", 24", 26" Transparent Spokeguards	£120
BRAKES*			
DCA7200		Pull Brakes	◆
DCA7210		Standard Push to Lock Brakes, Lever Bent	○
DCA7220		Standard Push to Lock Brakes, Lever Straight	○
DCA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre	(NC 7300, 7810, 7820, 7830) £74
DCA7240		Active/Scissor Brakes	(NC 7300) £74
DCA7250		Light/Compact Active, Aluminium Brakes	(NC 7300) £129
BRAKE LEVER EXTENSION			
DCA7300		Brake Lever Extension	£47
SAFETY OPTIONS			
DCA7500		Anti-tipper Right Maximum of 2 adoptions per side, includes Anti tips, Tipper aid and Cane holder	£116
DCA7510		Anti-tipper Left Maximum of 2 adoptions per side, includes Anti tips, Tipper aid and Cane holder	£116
DCA7520		Anti-tipper Both sides Maximum of 2 adoptions per side, includes Anti tips, Tipper aid and Cane holder	£232
DCA7600		Transportation Kit Includes Standard Pelvic Belt and Tiedown Labels	(NC 7800) £70
DCA7800		Standard Pelvic belt	£70
BODYPOINT BELTS			
DCA7810		BODYPOINT POSTURE BELT EVOFLEX SMALL	£77
DCA7820		BODYPOINT POSTURE BELT EVOFLEX MEDIUM	£77
DCA7830		BODYPOINT POSTURE BELT EVOFLEX LARGE	£77
OTHER OPTIONS			
DCA8000		Rear frame protection, plastic guard attached to underside of frame	£18
DCA8010		Front frame protection, 2 piece removable with hook and loop	£35
DCA8100		Tipper-aid, Right	(NC 8200) £44
DCA8110		Tipper-aid, Left	(NC 8210) £44
DCA8120		Tipper-aid, Both sides	(NC 8200, 8210) £87
DCA8200		Cane holder, Right	£62
DCA8210		Cane holder, Left	£62
DCA8450		TABLE TRAY	£319
DCA8500		Air pump	£12
DCA8550		Tool kit	◆
DCA8680		24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included	£304
BODYPOINT ACCESSORIES			
DCA8800		BODYPOINT CALF STRAP	£54
DCA8850		BODYPOINT MOBILITY BAG	£52

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DCA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable
DCA1300	280 mm Single cross brace
DCA1305	300 mm Single cross brace
DCA1700	370 mm
DCA1705	380 mm
DCA1710	390 mm
DCA1715	400 mm
DCA1720	410 mm
DCA1725	420 mm
DCA1730	430 mm
DCA1805	370 mm
DCA1810	380 mm
DCA1815	390 mm
DCA1820	400 mm
DCA1825	410 mm
DCA1830	420 mm
DCA1835	430 mm
DCA1840	440 mm
DCA1845	450 mm
DCA1850	460 mm
DCA1855	470 mm
DCA1860	480 mm
DCA1865	490 mm
DCA1870	500 mm
DCA1900	200 mm
DCA1902	210 mm
DCA1932	360 mm
DCA1934	370 mm
DCA1936	380 mm
DCA1938	390 mm
DCA1940	400 mm
DCA1942	410 mm
DCA1944	420 mm
DCA1946	430 mm
DCA1948	440 mm
DCA1950	450 mm
DCA1952	460 mm
DCA1954	470 mm
DCA1956	480 mm
DCA1958	490 mm
DCA1960	500 mm
DCA1962	510 mm
DCA2140	Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below)

DCA2155	Angle adjustable backrest, adjusted 82°
DCA2160	Angle adjustable backrest, adjusted 86°
DCA2165	Angle adjustable backrest, adjusted 90°
DCA2170	Angle adjustable backrest, adjusted 94°
DCA2175	Angle adjustable backrest, adjusted 98°
DCA2180	Angle adjustable backrest, adjusted 102°
DCA0000	Without angle adjustment. select with options 2110 fixed and 2120 lumbal backrests, supplied with fixed angle
DCA2210	Lightweight back upholstery, no tension adjustment
DCA2400	Without Push Handles
DCA2430	Push handles height adjustable, rear set (maximum height of push handles is rear seat height plus 600mm)
DCA2490	Headrest
DCA2595	Backrest Stabilization Bar, Swingaway and Removable
DCA2710	Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk
DCA2720	Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk
DCA2748	Contour Visco NG Cushion 100mm High - High Risk
DCA2750	Vicair Active O2 lightweight cushion 60mm High - High risk
DCA2810	Fluid supplement sac (Pair)
DCA2820	Tapered Wedge Pair
DCA2830	Abductor
DCA2840	Lateral Side Wedge
DCA2850	Anterior Front Wedge
DCA3000	Legrest 70° RH
DCA3010	Legrest 80° RH
DCA3020	Amputee legrest RH
DCA3030	Angle adjustable/Elevating legrest - RH, including angle, width and depth adjustable plastic footplate
DCA3040	Without Legrests RH, Select with fixed front frames
DCA3100	Legrest 70° LH
DCA3110	Legrest 80° LH
DCA3120	Amputee legrest LH
DCA3130	Angle adjustable/Elevating legrest - LH, including angle, width and depth adjustable plastic footplate
DCA3140	Without Legrests LH, Select with fixed front frames
DCA3260	2 pc. Plastic flip up
DCA3280	1 pc. Carbon flip-up to the left, angle and depth adjustable
DCA9999	No footrest (Select when ordering Amputee brackets both sides)
DCA3300	Calf strap mounted on the frame
DCA3305	Calf strap mounted on bracket
DCA3310	Heel loops
DCA3430	Carbon mudguards, fixed
DCA3500	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH
DCA3510	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH
DCA3520	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included
DCA3530	Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH
DCA3540	Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH
DCA3600	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH

DCA3610	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH
DCA3620	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH
DCA3630	Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH
DCA3640	Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH
DCA4190	Bold Decal (Full coloured transfer on frame)
DCA4195	Tag (Small Küschall brand transfer on front of frame)
DCA4200	Black
DCA4210	Dark Grey
DCA4220	Light Grey
DCA4230	Red
DCA4240	Blue
DCA4250	Gold
DCA1111	No bold colour, select only with Option DCA4195 Tag or DCA4197 Without Tag
DCA4300	Black
DCA4310	Red
DCA2222	No tag colour, select only with Option DCA4190 Bold Decal or DCA4197 Without Tag
DCA5550	150mm x 30mm (6" x 1.2") Comfort Black hub and tyre
DCA5650	175mm x 30mm (7" x 1.2") Comfort Black hub and tyre
DCA6000	One Arm Drive system Right side (with light wheel 7 bar tyres, alu handrims)
DCA6010	One Arm Drive system Left side (with light wheel 7 bar tyres, alu handrims)
DCA6040	Rear wheel extension, 6 rear wheel positions
DCA6100	1st Position
DCA6130	4th Position
DCA6140	5th Position
DCA6150	6th Position
DCA6300	Quad quick release/tetra clip
DCA6400	22"
DCA6410	24"
DCA6420	25"
DCA6430	26"
DCA6510	22", 24", 26" Küschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim
DCA6520	24", 25" High Performance Wheel Black hub and rim
DCA6550	24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black
DCA6580	22", 24", 26" Wheels with Attendant Drum Brake
DCA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DCA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
DCA6660	24" x 1 3/8" Solid Profile Black Tyre
DCA4444	Select with options 6000 and 6010 One Arm Drive wheels, supplied with specific tyre
DCA6700	22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DCA6710	24", 26" Steel, short tab with 1 hole, mounted at 3cm distance
DCA6720	24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DCA6740	24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip
DCA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands

DCA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion
DCA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim
DCA6780	22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DCA6790	22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting
DCA6840	Camber 3°
DCA6900	22", 24", 26" Transparent Spokeguards
DCA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DCA7300	Brake Lever Extension
DCA7600	Transportation Kit Includes Standard Pelvic Belt and Tiedown Labels
DCA7800	Standard Pelvic belt
DCA7810	BODYPOINT POSTURE BELT EVOFLEX SMALL
DCA7820	BODYPOINT POSTURE BELT EVOFLEX MEDIUM
DCA7830	BODYPOINT POSTURE BELT EVOFLEX LARGE
DCA8200	Cane holder, Right
DCA8210	Cane holder, Left
DCA8450	TABLE TRAY
DCA8680	24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):